

HEAT THERAPY

COMBINED HEAT THERAPY UNITS BTL-25

The BTL-25 L and BTL-25 XL bring a smart and effective solution to heat therapy. The devices combine two of the most popular kinds of heating therapy - paraffin baths and hot packs heating. Switching between therapies is easy by simple addition or removal of the stainless steel bath. The units are controlled by the digital thermo-regulator which provides exact temperature level and guarantees service effectiveness.



BTL-25 L Combined paraffin baths and hot packs heater



BTL-25 XL Combined paraffin baths and hot packs heater

Technical parameters	BTL-25 L	BTL-25 XL
Ordering No.	P625.021v100	P625.023v100
Regulator 30-90°C	digital	digital
Max. output	1500W (700W)	2×1500W (2×700W)
Water outlet	yes	yes
Voltage	230V	230V
Capacity - paraffin	20l	2×20l
Capacity - hot packs	5 hot packs (24l)	10 hot packs (2×24l)
Dimensions	330 × 540 × 290mm	660 × 540 × 290mm
Paraffin melt time	approx. 1 hour	approx. 1 hour

